



CLAIR DE LUNE  
DENTAL LABORATORIES

knowledge. precision. æsthetics. **function.**

# The Catalog



A Face by Moonlight #8.

# Our philosophy

“As long as the sun shines, one does not ask for the moon,” says an ancient Russian proverb. For a patient—a person—who is coping with the loss of a natural tooth and the adjustment to a prosthesis, the truth of the proverb stings the mouth and the soul.

We find satisfaction in creating a restoration of such beauty and natural grace that its new owner forgets... and smiles again.



A Face by Moonlight #9.

# Removable

Both denture packages feature a highly natural tooth arrangement, as well as an option for custom stained, shaded, and textured gingival bases. Our cast reinforcing framework is made from cobalt-chromium alloy for both strength and lightness.

A processed baseplate ensures accurate fit and is fabricated from high-impact denture resin. This procedure is highly recommended.

Another highly recommended option is an “Avoid Embarrassment” denture (AED) to give your patients a peace-of-mind backup during relines or repairs.

## **THE NEUTRAL ZONE DENTURE**

Our premier denture includes gothic arch tracers, an æsthetic blueprint, a neutral zone rim, one wax try-in setup, acrylic or ceramic teeth, processing, and finishing. Options include processed baseplates, glass-fiber reinforcement, cast reinforcement, stained and textured gingiva, gold posterior occlusals, and an AED.

## **THE DELUXE DENTURE**

Includes an æsthetic blueprint, one wax try-in setup, acrylic or ceramic teeth, processing, and finishing. Options include processed baseplates, glass-fiber reinforcement, cast reinforcement, stained and textured gingiva, gold posterior occlusals, and an AED.



A Face by Moonlight #11.

# Removable

## **THE CAST PARTIAL DENTURE**

Includes case design, a cobalt-chromium framework, one wax try-in setup, acrylic teeth, processing, and finishing. Options include stained and textured gingiva, and posterior gold occlusals.

## **TREATMENT DENTURE, TREATMENT PARTIAL DENTURE**

A treatment denture is a bandage for immediate post-extraction delivery and includes economy teeth and surgical guides. Our treatment partial includes wrought wire clasps. Options include glass-fiber reinforcement and cast clasps.

For those doctors who follow the branching technique, our treatment denture is ideal for temporization.



A Face by Moonlight #10.

# Fixed

We use only certified high-noble alloys and Shofu Vintage Halo porcelains and work under 20x magnification.

## **PFM CROWN, BRIDGE**

Our goal is simple: the patient cannot tell the difference between the natural and the prosthetic tooth in æsthetics and function. Typically, we create PFMs in conjunction with hybrid fixed-removable cases.

## **GOLD CROWN**

We create functional restorations with realistic morphology. For those doctors who follow the Tucker techniques, we use JRVT gold alloy with smoothed morphology. Typically, we create gold crowns or bridges in conjunction with hybrid fixed-removable cases. We also offer gold posterior occlusals with our dentures.



Two Faces by Moonlight #1.

# Combination

Clair de Lune specializes in hybrid fixed bridge-removable partial combination cases. For the best results, we work with you to plan the case before tooth preparation or implant placement. Options include acrylic, ceramic, or gold teeth.

We excel at implant and natural root-form locator attachment overdentures, hybrid fixed-detachable dentures, dentures with telescoping abutments and secondary frames, as well as bar-retained implant overdentures. Many other attachment systems are available.



A Face by Moonlight #12.

# Training

We offer free consultation and case design for everything we make. In addition, we offer face-to-face training for many procedures.

## **RELINES AND REPAIRS**

Learn how to offer same-day walk-in service.

## **CUSTOM TRAY FABRICATION**

Learn how to make your own trays.

## **TEMPORARY CROWNS**

Far beyond basic temporaries, discover guided tissue management as well as functional and æsthetic benefits of “Stage 1” indexing and “Stage 2” provisionals.

## **FACEBOW MOUNTING**

Learn how to properly mount a facebow to increase the accuracy of your articulated models and the resulting prostheses.

## **ANTERIOR TOOTH SETTING**

Learn to make chair-side æsthetic positioning of denture teeth during the bite registration or wax try-in phase.



Two Faces by Moonlight #2.

# Our policy

If you don't see a suitable package or service, please ask. We are happy to tailor any package or create a new one to meet your needs.

## **PRICES**

Please inquire about our prices. Alloys, teeth, and attachments are charged at market prices.

## **SHIPPING**

Shipping is extra.

## **WORKING TIMES**

We schedule all cases based on their complexity and our workload. You are always welcome to pre-schedule cases.

In general, we need 5 business days for each stage of removable denture cases, 5 days for each stage of fixed work, and 10 days for a cast removable partial denture framework.

## **REMAKES**

Inferior work does not leave our lab. We choose to work with doctors who share the same commitment to excellence. When a mistake happens, we use the opportunity for discussion and mutual learning. We offer a 50% discount for the remake.

## **PAYMENT**

We accept checks and all major credit cards. We ask for a credit card on file and submit invoices payable by check or credit card.

At Clair de Lune Dental Laboratories, we produce beautiful and accurate removable, fixed, and combination dental prosthetics with the finest bio-compatible materials because we care about our community, our work, and your patients' health.

**[www.cdldl.com](http://www.cdldl.com)**